

Issue Date:

Review Before:

Title: SOP FOR STORAGE OF STERILE ARTICLES IN STERILE STORE.

OBJECTIVES:

To lay down a procedure for the storage of sterile articles in sterile store.

SCOPE:




This procedure is applicable to store the articles after sterilization process in sterile store.

RESPONSIBILITY:

Assistant Supervisor, Sterilisation Technician, CSSD Attendants.

PROCEDURE:

1. After sterilization process the sterilized articles are checked whether the sterilization process is OK or not by checking the labels of articles.
2. After sterilization process shutter of drum is closed immediately after taking out from the sterilizer.
3. Overshoes, cap, mask and proper dress should be wear before entering in sterile store.
4. Sterilized articles are stored in sterile store at their respective places on slabs & racks.
5. Sterile articles are issued from issuing window attached to sterile store.
6. Fumigation should be done time to time.

	Prepared By	Checked & Approved by	Authorized By
Designation	St. technician	Asstt. Supervisor	Officer In charge
Name	Anil Kumar	Yogesh Kumar	P.K. Verma
Sign/ Date			 7/10/11